诚实守信 合作共赢
Honesty Trustworthy, win-win Cooperation

PENG RUI ENTERPRISE

BIOLOGICAL/PHARMACEUTICAL MACHINERY
EQUIPMENT/SPECIAL ALBUM
生物制药机械设备 专集

› 符合FDA标准/满足GMP要求
Accord FDA Standard/meet GMP Requirement

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
提取车间 Extraction workshop

饮料车间 Beverage workshop

正锥型提取罐
Normal taper type extracting tank

直筒型提取罐
Straight tube type extracting tank

蘑菇型提取罐
Mushroom type extracting tank

倒锥型提取罐
Inverted cone type extracting tank

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This machine is suitable for pharmacy, biology, food, chemical industry, often used in Chinese medicine effective ingredients and balmy oil extraction and concentration of traditional Chinese medicines, is the common equipment.

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
A. Strictly according to pressure vessel code, precision polishing, make internal and external comply with GMP. B. Equipment is mainly by the cans, away PaoQi, condenser, cooler, the oil-water separator, filters etc. C. The main form of the equipment has straight tube type, is cone type, grinding, she reversed cone of such as type, can according to need equipment mixing device. D. Dregs door have Duo Zhong form for the user choice, there is double cylinder type, many cylinder type, spiral, tooth jewels. E. According to customer needs can increase when the liquid filtration area to improve production efficiency.

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RH HOT REFLUX EXTRACTION ENRICHMENT UNITS (LOW TEMPERATURE)

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

This equipment is suitable for extracting effective components of Chinese herbal medicine and many plants, concentration, solvent recovery and collection of aroma.

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This equipment can be used for dynamic continuous ultrasonic extraction. At the same time, it can be equipped with a large capacity stainless steel rotary evaporator specially developed for schools and scientific research institutions to concentrate the medicine extracted by the ultrasonic extraction machine. The equipment can directly concentrate the liquid medicine to the state of extract. This set of equipment is very suitable for universities, research institutes, enterprises and institutions, laboratory R & D departments and pharmaceutical factories for multi variety and small batch production. Ultrasonic extraction (also known as extraction) with its unique advantages of low extraction temperature, high extraction rate and short extraction time has been applied to the extraction of effective content of Chinese medicinal materials and various animals and plants by innovative people. It is a modern high-tech means to replace the traditional shearing process to achieve efficient, energy-saving and environmental protection extraction.

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1. High extraction efficiency: the unique physical characteristics of ultrasound can promote the cell wall breaking or deformation of plant tissue, so that the effective components of traditional Chinese medicine can be extracted more fully, and the extraction rate is significantly increased by 50-500% compared with the traditional technology;

2. Short extraction time: the best extraction rate can be obtained in 24-40 minutes by ultrasonic wave. The extraction time is greatly shortened by more than 2/3 compared with the traditional method, and the processing capacity of raw materials is large;

3. Low extraction temperature: the optimal temperature for ultrasonic extraction of Chinese herbal medicines is 40-60 °C, which can protect the active ingredients in Chinese herbal medicines which are unstable in heat, easy to hydrolyze or oxidize, and save energy and energy greatly;

4. Wide adaptability: ultrasonic extraction of Chinese medicinal materials is not limited by the polarity and molecular weight of components, and is suitable for the extraction of most kinds of Chinese herbal medicines and various components;

5. There are few impurities in the extracted liquid, and the effective components are easy to be separated and purified;

6. The operation cost of extraction process is low and the comprehensive economic benefit is remarkable;

7. The operation is simple and easy, and the equipment maintenance is convenient.

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The results showed that the extraction rate was increased by 50% - 500%; the extraction time (min) was shortened by more than 2/3; the extraction temperature was 40-60 °C; the National Research Center for modern engineering and technology of traditional Chinese medicine has made a variety of medicinal materials, compound compounds and various chemicals, such as flavonoids, etc. The ultrasonic extraction of saponins, terpenoid quinones, coumarins, lignans and alkaloids was compared with the traditional extraction methods. The relationship between the extraction rate and ultrasonic time, power and frequency, and the relationship between solvent concentration and extraction rate were determined. The optimal process conditions were established and valuable first-hand information was accumulated. Thcg series ultrasonic extraction tank is a new ultrasonic extraction technology based on the original extraction tank. It can be divided into two types: internal and external. In animal and plant extraction, ultrasonic cycle extraction technology has been widely used. Ultrasonic extraction is to use the cavitation effect of ultrasound to increase the penetration of solvent, improve the dissolution rate and dissolution times of drugs, so as to increase the diffusion of material components, shorten the extraction time and accelerate the extraction process.

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This equipment is suitable for extracting effective components of Chinese herbal medicine and many plants, concentration, which can realize the solvent recovery and aromatic oil collection. Especially suitable for hospitals, pharmaceutical factories, scientific research units to develop products and small batch production.

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Equipment descriptions

The single-effect external circulation vacuum concentrator is composed of heater, evaporator, defoamer, condenser, cooler and liquid receiving tank. The parts that contact materials are all made of stainless steel. The equipment is applicable to evaporation and concentration of liquid materials in pharmacy, food, chemical industry and light industry etc. industries. The concentrating time is short and the evaporating speed is high, which can better protect the heat sensitive materials from being damaged.

WZH 单效外循环真空浓缩回收器
WZH single effect external circulation
Vacuum recovery concentrator

WZ 单效外循环真空浓缩器
WZ single effect external circulation
vacuum concentrator

WZS 单效外循环真空浓缩收膏器(自动防爆型)
WZ single effect external circulation
vacuum cream reclaim concentrator

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The series of single-effect evaporators have three forms

I WZ single effect external circulation vacuum concentrator is the ordinary type, which is used for evaporation of water solution.

II WZH single-effect external circulation vacuum recovery concentrator is the alcohol recovery type, which is used for evaporation of alcohol solution. It has higher cooling area to ensure the alcohol recovery. Its recovery strength reaches 70~85%, the alcohol consumption rate is less than 5% and the alcohol temperature after the recovery is less than 40C.

III WZS single-effect external circulation vacuum ointment collection concentrator is used for single-effect evaporation and thick ointment collection. After the normal single-effect evaporation, the valves connecting the spraying tube and circulating tube in the heating chamber and evaporating chamber will be closed. The liquid materials will fully enter the evaporating chamber and through the jacket heating in the evaporating chamber, normal pressure type or pressure relief type concentration will be carried out. The concentration specific gravity of materials will be increased to 1.35-1.45.

注：根据客户要求推出防爆型全自动型浓缩器，全程温度控制、真空自动调节、液位控制及自动排水等并能时时记录保存数据
Note: The explosion-proof type automatic concentrator is developed according to customer requirements. The whole-course temperature control, vacuum automatic adjustment, liquid level control and automatic water discharging can be realized and data can be recorded and saved in a real-time manner.

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Application range: applicable to low temperature vacuum concentration of heat sensitive materials, especially applicable to concentration of traditional Chinese medicine, western medicine, starch sugar and other water solution.

Working principle: After being compressed through hot pump, the first-effect secondary steam is supplied to the first-effect heating, so that large amount of steam is saved and production cost is reduced.

Equipment composition: first-effect heater, first- effect evaporation chamber, hot pump, second effect heater, second-effect evaporation chamber, condenser, liquid receiving tank, a set of pipe and valve in the system.

Temperature index: live steam pressure 0.5~0.7MPa; first-effect heating temperature 90°C; first-effect evaporation temperature 75°C; second-effect heating temperature 75°C; second-effect evaporation temperature 75°C; second-effect evaporation temperature 55°C.

Steam consumption: The steam-water ratio is 0.33-0.37kg steam/kg water. It can save 34% energy and 30% water than the dual-effect concentrator without hot pump.

技术参数 Technology parameter

项目 单位		带热泵中药双效浓缩装置 Tradotomal Chinese medicine dual-effect concentrating device provided with hot pump	带热泵中药双效浓缩装置 Ordinary traditional chinese medicine dual-effect concentrating device
高压引射蒸汽压力 High pressure ejection steam pressure	Mpa	0.5~0.6	0.1~0.25
热能压缩器噪声 Noise of hot energy compressor	dB(A)	80~85	
理论汽水比 Theoretical steam-water ratio	kg汽/kg水	0.30	0.50
实际运行汽水比 Actua steam-water ratio	kg汽/kg水	0.33~0.35	0.65~0.75
第一次	壳程真空度 Shell side vacuum degree	MPa	0.040~0.030
	加热蒸汽温度 Heating steam temperature	°C	85~90
	蒸发室真空度 Evaporation chamber vacuum degree	MPa	0.068~0.060
	沸点 Boiling point	°C	70~75
第二次	壳程真空度 Shell side vacuum degree	MPa	0.068~0.060
	加热蒸汽温度 Heating steam temperature	°C	70~75
	蒸发室真空度 Evaporation chamber vacuum degree	MPa	0.089~0.082
	沸点 Boiling point	°C	52~57
各效加热室凝结水排出方式 Discharging method of condensed water in chambers of various effects		等压排水串级自蒸发, 可回收凝结水部分热量。With isopistic pressure water discharging serial self-evaporation, some heat of the condensed water can be recycled.	
冷却水量 Cooling water volume	进入温度25出水温度35 Water inlet temperature 25 water outlet temperature 35 °C	取此装置冷却水量为1 The cooling water volume of the device is taken as 1.	相同蒸发量, 冷却水量为1.43 Under same evaporation, the cooling water volume is 1.43.
操作性 Operability		操作简便, 调整及观察容易 Simple operation and adjustment and easy observation	操作简便, 调整及观察容易 Simple operation and adjustment and easy observation
经济性 Economical efficiency		占地面小, 设备造价较高, 使用成本低, 适合热敏性物料 Small land area, higher equipment cost and low using cost Suitable for heat sensitive materials	占地面小, 设备造价低, 使用成本高 Small land area, lower equipment cost and high using cost
增加设备投资成本 Increasing equipment investment cost		1.50	以此项成本为1 Taking the cost of the item as 1

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This equipment is suitable for pharmaceutical, foodstuff, chemical and other industries of material liquid concentration, distillation and recycling the organic solvent, etc. As a result of reduced pressure concentration, and concentration time is short, for heat-sensitive materials won't destroy the effective ingredients.

This equipment is mainly composed of concentrator subject, condenser, liquid separator, the liquid barrel of four parts, the equipment and material contact parts made in 0Cr18Ni9 stainless steel, and has good corrosion resistance, durability, and conform to the requirements of the drug and health.

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The device is suitable for alcohol recovery in pharmacy, food, light industry, chemical industry and etc., with the aim to reduce production cost and enhance economic benefit. The device can distill the thin alcohol of 30-50% to 90-95%. If the requirements to the finished alcohol strength is high, the reflux ratio can be made higher, but the output of finished products will be correspondingly reduced. The alcohol content of discharged residual liquid is less than or equal to 0.1%. It is a device specially for alcohol recovery.

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This equipment by column reactor, minus its, condenser, cooler, buffer tank high storage tank, six parts, parts contacted with material using high quality stainless steel SUS304 or SUS316L manufacture, has good corrosion resistance, inferior smooth processing, on the surface of the GMP standard, this device can will be 30 degrees~50 degrees thin alcohol distillation of to 93 degrees to 95 degrees, residual liquid containing alcohol degrees, low emissions and energy saving, environmental protection, to reduce production cost and improve the efficiency of the advantages.

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The equipment is suitable for alcohol sediment in Chinese medicine, oral liquid, foodstuff, health-protective products and chemical industry. Alcohol sediment is developed from traditional technology, mainly used for freeze Chinese medicine liquid being fried and condensed or to operate warm alcohol sediment, also can be used to operate water sediment after Chinese herbs being extracted.

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Equipment characteristics

The device includes the ellipse seal head with clip sleeve, tapered base of cylinder, inner-installed three paddle miller and special finely regulated revolve-out liquid pipe, etc. The mode for water-decocting and concentration of China traditional drug as follow. First fill the alcohol to form mixture, secondly under low temperature freeze to make deposition to separate the solid from the liquid to improve the purity and cleanness of the extraction liquid of Chinese traditional drug. The device is made of stainless steel of SUS304. The high pressure automatic spray & cleaning system is available. It is available to fill the low temperature of cooling water through the clip sleeve. The mechanical seal in the transmission part and the explosion-proof motor can ensure the security of the production process, meeting with the GMP standard. The surface is treated with surface mirror or darkness treatment.

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Tanlet vessels, creating aseptic status through specialization

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Constantly improving production details of Tanlet

1、After the order for our product is made, the steel plate purchased does not contact with the cement ground. The wood board or rubber cushion is adopted for support during all working procedures including cutting, polishing and forming. In this way, the steel plate surface is protected from any damage. Otherwise, once it is scratched, the polishing work will be added and more polishing work will reduce the thickness of the steel plate.


2、During the welding for product forming, the plasma automatic welding machine imported from Germany is adopted. The intercrystalline granules after the welding are more even with higher strength as compared to the ordinary argon arc welding machine.

3、The interior surface of the equipment achieves a roughness of Ra<0.4 μm. The round corner is hydraulically shaped through adopting Japanese style flat plate to ensure there is no dead corner. The polishing technology adopted is automatic polishing. We have independently developed internal polishing machine and external polishing machine with advanced technology. For some devices with higher requirements, it is required that the welding seam position should not be found through naked eyes or hand feel after the polishing. The exterior surface adopts upper and lower seal heads with mirror polishing to avoid light contamination. The middle tank adopts matte treatment or is manufactured by abrasive plate.

4、The product is carefully inspected, even the detail is paid much attention to. For example, the pipe end height is required to the nearest mm and the verticality of pipe end is required to be within 2 degrees.

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不锈钢容器系列 Stainless steel tank series

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This equipment is widely used in the dairy industry, food, pharmaceutical, chemical, beverage and other industrial sectors, hot and cold cylinder as the heating, cooling, heat preservation, sterilization or slurry storage is a necessary equipment, especially suitable for cold and hot cylinder of pasteurized milk and heat preservation, also is the necessary equipment in ice cream production line.

Volume: 50L, 100L, 200L, 300L, 500L, 600L, 1000L ~ 1000L, stirring speed, 15~ 120 R/min speed.

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This equipment is widely used in food, beverage, biotechnology, pharmaceutical and other industries of material using enzymes as catalyst, to break down all kinds of organic polymer, make it into a small organic molecules, facilitate the organism assimilation. Tanks with interlining, can pass into the steam or cooling water respectively for indirect heating or cooling, make the material under the action of enzyme, using blender running at optimum temperature and pH conditions fully enzymolysis reaction. Volume: 50L, 100L, 200L, 300L, 500L, 600L, 1000L ~ 1000L, stirring speed, 15~ 120 R/min speed.

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This ingredients tank also known as ingredients cylinder, cylinder, mixing barrel, mixing tank, energy saving, durable corrosion, strong production capacity, easy clean, simple structure, etc; Designed for food (milk) products, beverage, pharmaceutical, bio-engineering, wine, spices, flavoring, daily chemical industry design and manufacture of stirring and mixing equipment. Tanks with interlining, can pass into the steam or cooling water, respectively, for heating or cooling, so more than two kinds of material (such as powder or additives + solution, etc.) with the help of agitator operation under the condition of optimum temperature for fully mixing, stirring. Volume: 50L, 100L, 200L, 300L, 500L, 600L, 1000L ~ 1000L, stirring speed, 15~ 120 R/min speed.

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When products need agitation, according to your actual process requirements Tanlet company makes design and correctly choose blade type which can result in agitation results of optimization, so that it lets products quicker and equaler movement.

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
It belongs to axes flow with high speed circulation and little cut force. It uses baffles or draft tube to make axes more strength circulation, applies on mixture, transfer-heatAnd solid suspend Low-viscosity soluble liquids. Sanitary grade equipments must reduce dead Corners and be easy to clean, So it does not suit for adding baffles and draft tube. Yuanyue designs off-centre installation that solve axesflow and intense turbine and rotation agitation causes.

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There is all horizontal circle flow in different height of blade. it uses oblique or angle steel type blades that increases turbine around blades. Anchor type and frame type operates in laminar flow for increasing scope of agitation, according to requirements, adds vertical blades and beam which suits formixture of high-viscosity liquid or liquid including high solid.

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The crystallizer uses imported 304 or 316L to be fabrication. With inner wall of mechanical polishing or electrolytic mirror polishing, outer wall of all welded structural insulation. outer surface of mirror or matt.

The respirator of 0.22um is installed on the crystallizer.

The crystallizer can bear 121C high temperature sterilization and disinfection with sanitary pressure meter.

Outer connection of the crystallizer adopts sanitary quick connector.

The crystallizer has temperature monitor device, two sight lights, cleaning ball and inlet of Ni, temperature connection, etc. Seal of shaft adopts special

sanitary mechanical seal design which insures they are not polluted.

The crystallizer uses adjustment device of frequency conversion with big range adjustment of shaft speed.

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SS reactor tank is one of normal use reactor equipments in pharmaceutical and chemical fields. at certain temperature and pressure, through agitator mixes two kinds of volumes certain (or many kinds of) liquids, solids or liquids, promotes its chemical reaction, usually with thermal heat effect, Inputs required heat energy by heat exchanging device or shifts forming heat energy. Type of agitation more uses anchor or frame type which ensures equal and even mixture in short time.

Speed of agitation more matches frequency conversion with temperature of display and control

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REACTING KETTLE AND BLENDING TANK IN FOOD FIELD

专用美拉得反应釜Special Meaction Type Reaction Kettle

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It is the reaction kettle for such substances as protein, fat, amino acid, peptide, sugar and etc. to produce the substances full of strong and intense fragrance of cooked meat after thermal reaction under high temperature. It is a special device for producing substances full of fragrance of cooked meat. It can be used for other reaction or used as a small preparation tank. The equipment has its special blender, and is provided with condenser for heating and reflux. The jacket can be connected with steam for heating, and with cooling water for quick cooling. Heat preservation layer is provided outside.

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
This machine is suitable for sugar, starch kind of protein powder, milk powder, kind of product and emulsifier, stabilizer, thickening agent to mix, usually in the formula milk powder, yogurt, all kinds of beverages, ice cream and other ingredients in the process.

By emulsion can block, the stator and rotor, mechanical seal, emulsification head, motor, rack and other components, the equipment is the most advanced domestic currently the emulsifying equipment, reasonable structure, convenient cleaning, convenient operation, good appearance.

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Main structure and configuration of high shear emulsion tank

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The shear head adopts claw type and two-way suction structure, which avoids the dead corner and vortex phenomenon caused by the upper material difficult to inhale. The high-speed rotating rotor produces strong shear force, which makes the shear rate higher and the shear force stronger. Under the action of centrifugal force generated by the rotor, the material is thrown into the narrow and precise gap between the stator and rotor from the radial direction. At the same time, it is subjected to centrifugal extrusion, impact and other forces, so that the materials are fully dispersed, mixed and emulsified.

Note: if it is used in vacuum or pressure vessel, the motor support of shearing machine shall be equipped with corresponding mechanical seal device.

高剪切乳化罐
HIGH SHEAR EMULSIFICATION TANK

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Widely used in big infusion and injection, solution, biological products, oral liquid, fermented liquid materials such as the preparation

Volume: 10L-300L, this equipment structure design is beautiful, convenient operation, supporting feet adopt universal wheel design, health, easy clean, flexible mobile firm; Inside, outside surface by the health level mirror polishing treatment in accordance with GMP standard of quality.use of mechanical seal, explosion-proof electric motor to ensure production process safetyGMP standard of quality.

实验室反应釜
LABORATORY REACTION KETTLE

电加热搅拌罐
ELECTRICAL HEATING MIXING TANK

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Suitable for: pasteurized milk, yogurt, fancy milk, juice, tea beverages, beer, ice cream, milk beverage, concentrated juice, wine etc.

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Pasteurising technology parameters: 4°C to 75 C and homogeneous 85~95 C (15S)-15 C-4 C

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

The kettle, kettle by kettle body cover, the clip set, mixer, transmission shaft sealing device, device, supporting etc,Device equipped with CIP cleaning rotating spray cleaning ball head, the WenKong, pressure taps, explosion-proof depending on the light, as a mirror, quick release open mouth and so on, ensure convenient operation, use reliable. The pollution, the noise of the production environment characteristics, such as small this equipment all use SUS304/SUS3 16L stainless steel production, both internal and external surface sanitation mirror polishing or drawing processing, in accordance with The transmission parts

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This equipment is used in medicine (material workshop, synthetic workshop), biological engineering, fine chemical, chemical synthesis, food industry etc of hydrolysis, and crystallization, distillation, evaporation, heating mixed.

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THE OPENING TYPE PRESSES THE FILTER QUICKLY STAINLESS STEEL FRAME FILTER

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The Opening Type Presses The Filter Quickly:

Quick open connection, easy to clean; design to be pressurized, pump transfer to be filtered, and vacuum filtering; imported filtering mesh bag, high filtering precision, durable; machine connection mechanical prolonged, circle arc transition; polish degree of inner cell and outer shell of 0.4μm; and high precision pressure application, anti-vibration, water proof.

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DL Series Stainless Steel Frame Filter:

The machine is stainless steel multilayer frame filter. It is fit for enclosed filter in liquids with concentration 50%, low viscosity and little slag so as to meet precision filter and half-precision filter requirements like refining, sterilization and clarification. By directly using micro porous-membrane, material can achieve aseptic filtration without passing micro porous-membrane filter. The filter has big filter area, big flow, and wide application scope. So it is widely used in pharmacy, chemical engineering and food industries. While used for filtering injection liquids in pharmaceutical factories, it has particularly good effects.

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
Main Technical Parameters:

Filtration precision: 4 μm-50 μm
Filtration efficiency: 95 % -99.99%
Working pressure: ≤0.3MPa
Working temperature: ≤230C
Specification of Titanic stick: 中60X300mm

钛棒与流量对照表 Contrast form of Titanic stick and flow


型号Model	TL-1	TL-3	TL-5	TL-7	TL-9	TL-11	TL-19
规格芯Specification Core	1	3	5	7	9	11	19
流量FlowT/h	0.3	0.9	1.5	2	2.5	3	5
外形尺寸Dimension LxWxH(mm)	610x185	385x350x865	395x410x880	530x430x890	500x465x880	500x465x890	600x600x930

Equipment characteristics:

Titanium powder sintered filter core is adopted, usually used for prefiltration or interfiltration. This filter core has high precision, high mechanical intensity and is high temperature resistant and anti-corrosion, which is widely applied to medicine liquids decarbonization and gas filter. It can also replace sand core filter that is commonly used on large liquid infusion production lines. The shell of TL titanium rod filter of our company is made of SUS304 or SUS316L stainless steel. After used by customers for years, it has been proved to have good effects and technical indexes apply with GMP medical standard. The titanic stick is made of sintered titanium alloy powder. The products filtered by it have such characteristics as solid and liquid separation, particle removing, heat source removing, high clarity and good mechanical and physical performance. It is mainly used for decarbonization of medicine liquids, and it can also replace sand stick filter that is commonly used in large transfusion production. After used by customers for years, it has been proved to have good effects and all its technical indices comply with pharmacopoeia stipulations and GMP requirements.

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Plate and frame filter is a kind of filtering equipment, applied to various suspension of solid-liquid separation, wide applicable scope, applicable to the pharmaceutical, food, chemical industry, environmental protection, water treatment and other industrial fields.

Plate and frame filter is suitable for low concentration below 50% viscosity, less fluid slag content for airtight filtering to achieve purification, sterilization, clarify the purpose of the filter, etc.

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Diatomite filter is suitable for white wine, fruit wine, wine, health wine, beverage, syrup, soy sauce, vinegar, and medicine, biological, chemical and other clear filtration of liquid products. Filtration cycle length, high efficiency, stable turbidity, good sealing, compact structure, convenient operation, portable, easy to maintenance.

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The contact parts of Pengrui Health Pump series are SUS 304 or SUS 316L. Steel manufacturing, widely used in medicine, beverage, chemical industry, brewing and other industries, in line with the pharmaceutical industry GMP requirements, welcomed by pharmaceutical companies. SUS304 or SUS316L for pump body and impeller Made of dense casting or forging. The internal finish of the pump is polished to reach Ra 0.28 μm ~ 0.4 μm. Import and export adopt international standards ISO, DIN standards. Flow: 3T/h - 30T/h lift-out: 20m ~ 36m.

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


Pengrui people are well aware that the close relation between tube element, valve and cleaning equipment is important for the sanitation requirement of the whole set of technological facilities. Therefore, each part is elaborately processed, materials all adopt SUS304、SUS316L, internal surface has mirror polish with finish Ra≤0.28 μm, surface adopts mechanical polishing, sand grinding or surface polishing. Besides, International technical criteria such as ISO、DIN、IDF and SMS are adopted So the product series can go with GMP medical standards.

大输液配置工艺流程 Configuration process flow of infusion

小针剂配置工艺流程 Process flow of dispensing small injection

果醋工艺流程图
Process flow diagram